GeoSafari® STEREOSCOPE

EI-5303

8+

WARNING:

CHOKING HAZARD—Small parts.
Not for children under three (3) years.

Educational
Insights®

STEREOSCOPE DIAGRAM

1. Eyepiece (10x and 20x magnification)
2. Focus knob
3. Light
4. Stage
5. Light on/off button

WELCOME, YOUNG SCIENTIST!

Science is all about exploration and discovery. Your new microscope is an instrument with lenses for making small objects appear larger so they can be studied. There are millions of plants, animals, fossils, and rocks that you can scope out—in amazing detail—with this new tool.


In today's technological world, every field of science uses microscopes. Microscopes are used by biologists to study microorganisms, by geologists to study rocks and minerals, by archaeologists to study ancient artifacts, and by astronomers to study fallen meteorites.

The GeoSafari Stereoscope you just unpacked is perfect for viewing rocks and fossils. It is a special kind of microscope that allows you to look with both eyes at once to see fine detail on 3D objects. Get started with one of the included rocks and then keep exploring your world to find new objects to study!

BATTERY INSTALLATION INSTRUCTIONS


Note: Before installing batteries, check to make sure there are no eyepieces in the Stereoscope. If there are, remove them before turning the Stereoscope over to prevent them from breaking.

- Use a Phillips-head screwdriver to remove the battery door screws.
- Stereoscope requires 2 AA batteries, not included.
- Batteries must be inserted with the correct polarity.
- Do not mix old and new batteries.
- Do not recharge non-rechargeable batteries.
- Do not mix different types of batteries: alkaline, standard (carbon zinc), or rechargeable batteries.
- Only use batteries of the same or equivalent type.
- Remove exhausted batteries from the unit.
- Do not short-circuit the supply terminals.
- Only charge rechargeable batteries under adult supervision.
- Remove rechargeable batteries from the Stereoscope before charging.
- To prevent corrosion and possible damage to the product, it is recommended to remove the batteries from the units if they will not be used for more than two weeks.
- Please retain these instructions for future reference.

HOW TO USE YOUR GEOSAFARI STEREOSCOPE

1. Follow the battery installation instructions to install 2 AA batteries.
2. Choose an eyepiece. There are 10x and 20x eyepieces included. We recommend starting with the 10x eyepiece to find what you want to inspect and then increasing to the 20x eyepiece later to get a closer look.
3. Place the eyepiece in the microscope. The 10x eyepiece is installed when you first open your Stereoscope.
4. Press the on/off button to turn on the LED. This will light the object you place on the stage.
5. Place an object on the stage. We have included 12 rocks to get you started but you can inspect anything that fits on the stage!
6. Turn the focus knob to adjust the focus.
7. When you are done, press the button to turn off the LED.

ROCK SAMPLES

There are 3 main types, or classes, of rocks in the world and in this set. They are igneous, metamorphic, and sedimentary.

Igneous rocks are formed when magma (molten rock deep in the earth) cools and hardens. When magma escapes from inside the earth, it is called lava.

Metamorphic rocks are formed when one type or rock changes (metamorphoses) into a new type of rock due to extreme heat and pressure inside the earth.

Sedimentary rocks are formed when layers of material are added on top of existing rock over a long period of time. Sedimentary rocks are usually soft and may break apart or crumble easily.

*Rock colors and features may vary.

Obsidian is an igneous rock that forms when molten rock material such as lava cools very quickly. It usually has a very smooth texture and is also commonly referred to as volcanic glass. Obsidian is mostly black but can sometimes be brown, tan, or green and very rarely blue, red, orange, or yellow.

Obsidian


Granite is an igneous rock that is often mixed with magmatic rocks. Due to its strength and toughness, granite is often used for sculpting, construction, and design. Granite can be found in a variety of colors but is mainly white, pink and gray.

Granite

Pyrite is a mineral and has a bright, metallic finish. It is commonly called "fool's gold" since it looks so much like gold. Use your Stereoscope to find the gold flecks in the ore.

Pyrite Ore

Quartz is a widely found mineral that consists mainly of silica. It is the second-most abundant mineral in Earth's crust after feldspar. There are many varieties of quartz including gemstones like amethyst and rose quartz. The sand you see at the beach is even a type of quartz!

Milky Quartz

Calcareous Tufa

Tufa is a deposit of calcium carbonate, which is also what makes up an eggshell. It is a porous rock that often forms around mineral springs.

Muscovite

Muscovite is a mineral belonging to the Mica family. As a mineral, it can be found in all three types of rock. It can be sliced into thin, clear sheets and therefore has been used to make windows. It can also be crushed and used to add a glitter effect to paints or cosmetics.


Gneiss

Gneiss is a metamorphic rock that can be identified by its bands of varying composition. Gneiss often forms around convergent plate boundaries and forms under very high temperature and pressure.


Marble

Marble is a metamorphic rock composed most commonly of calcite. It is formed when limestone is subjected to heat and pressure inside the earth. Marble is commonly used for sculpting, construction, and design.

Graphite is a crystalline (organized) form of the element Carbon and it occurs naturally. It is often found in both metamorphic and igneous rocks around the world. It is a very soft, lightweight mineral and attaches to other materials very easily. It is also very heat-resistant, making it very versatile. Graphite is commonly used in pencils.

Graphite

Selenite

Selenite is a variety of gypsum and is characterized by a clear appearance. Gypsum is a soft mineral and is most often found in sedimentary deposits.

Scoria is an igneous rock with many round, bubble-like cavities. These cavities form when magma with lots of gas inside flows from a volcano. As the lava emerges, the gas begins to escape from the rock and the bubbles form. As a result, scoria is often lightweight. Scoria is usually black, gray, or sometimes even reddish brown.

Scoria

Limestone is a sedimentary rock and is composed mainly of calcium carbonate as the mineral calcite. Limestone usually forms in calm, shallow water from the accumulation of shell, coral, and algae. Because of this, limestone is often found in tropical locations with abundant reef life and activity. Chalk is a variety of limestone.

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FIND ALL OF THESE AND MORE AT EDUCATIONALINSIGHTS.COM

El-5305 GeoSafari® Omega Telescope

GeoSafari®

El-5301 GeoSafari® Micropro™
95-Piece Microscope Set

Developed in Southern California by Educational Insights.

All rights reserved. Made in China. ©Educational Insights, Gardena, CA, USA.

Learning Resources Ltd., Bergen Way, King's Lynn, Norfolk, PE30 2JG, UK.

educationalinsights.com

Please retain the package for future reference.

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⚠ WARNING: Cancer and reproductive harm. For more information go to www.P65Warnings.ca.gov.

**Educational
Insights®**